DIRECT SYSTEM

OCCLUSIONVD

Bio-Esthetic & Function


OCCLUSIONVD

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INNOVATION Thanks to their unique laser sintering process, the translucent enamel and occlusion veneers stand out due to their inorganic, ceramic-like surface vitrification and tempered composite body. The occlusions represent the anatomic foundations for single or complete reconstructions and for increasing the vertical dimension (OVD) in the posterior region. Hence, a functioning anterior cusped guidance can also be obtained using edelweiss „VENEERS“.

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INDICATIONS

Anterior and posterior restorations
Anatomical deformities
Erosion
Attrition
Lifting the vertical dimension (CMD)
Semi- direct and indirect restorations

ADVANTAGES

Minimally invasive
Clinically effective
Economical

The Progressive Technical Concept

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VENEER & OCCLUSIONVD the laser treated process combines the best of two worlds into one: a homogenous, inorganic and high gloss surface fused together with a thermally-tempered and dynamic composite core produces an optimal integration between function and esthetics. The difference is in it's similarity to nature.

smooth inorganic surface
without visible composite structure

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
SIZES Based on a study of all shapes and size variations of natural tooth anatomy, prefabricated and contourable universal occlusion shapes for the upper and lower arch were developed.

THE CHOICE of the tooth shape is performed with the available sizing guide. (OCCLUSIONVD - Sizing Guide). The sizing guide is positioned over the teeth to be restored, and the outline allows for proper selection of the best fitting occlusion.

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
MORPHOLOGY Ultra-thin, anatomical occlusions in three sizes (S, M & L). The glass component in the composite filler and laser treatment of the surface gives the material nature-like mechanical properties (biomechanics), and the material is also biocompatible (biology). Due to the natural morphology, the pre-fabricated Occlusions are very easy to incorporate into an existing occlusion.

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



FUNCTION VENEERS & OCCLUSIONVD allow minimally invasive biomechanical and economical treatment of the cause as well as protection of craniomandibular dysfunction (CMD). Furthermore malocclusion and deep bites can also get corrected minimally invasive. By covering the occlusal surfaces, usually in the mandible, the mandible is brought into the correct position in relation to the maxilla for the patient. This makes it possible to attain functioning guidance of anterior teeth and cuspids by using Edelweiss VENEERS.

LIFTING THE VERTICAL DIMENSION using Edelweiss occlusal surfaces. In most cases the mandible is appropriate. Can often be identified by the great difference in height from the cuspid to the premolar (33 to 34 and 43 to 44). However, to ensure that lifting, and protrusion if necessary, is equal on the left and right, this needs to be tested using situation models on the articulator! In the premolar region this can mean capping up to 5 mm. And up to a further 2 mm in the molar region.

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PROCEDURE The maxilla and the mandible models are articulated head-adequate in habitual teething. Then protrusion is performed as required and the supporting pin of the articulator raised. Partial bite splints with „stops“ are prepared for the mandibular quadrants. These are fixed temporarily in the mouth, in other words, provisionally. Already after a few days this will represent the ideal position of the mandible, possibly requiring reocclusion. On the basis of these bite splints the mandible is articulated anew versus the maxilla.

INDIRECT APPLICATION By lining the pre-fabricated occlusal surfaces on existing model occlusal surfaces with composite, one can simply create a new functioning biomechanical, balanced and smoothly working masticular apparatus without irritations. Then the Occlusions are transferred from the model to the patient's occlusions and attached via bonding and composite.

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DIRECT PROCEDURE for advanced professionals: Instead of lining the occlusions on the model in the articulator, a bite registry material (silicone) is used to make an occlusal impression which accurately portrays the clearance between maxilla and mandible in the articulator. Using the fixed silicone bite register, direct application of the prefabricated occlusion can be performed. The silicone bite register is in each case placed unilaterally (left or right) on the patient's occlusion. In other words, when the left quadrant is restored with occlusions, the right quadrant fixates the height and position of bite elevation and thus the final position of the composite-lined and bonded occlusions on the left when biting down (replaces the supporting pin when closing the articulator).

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Once the left quadrant has been restored, the right quadrant is treated. Now one can optionally prepare a further silicone register from the quadrant already restored with occlusions (left) in combination with the silicone register (right) already prepared with the articulator, which additionally guides and fixates the already attained occlusional relationship between the maxilla and the mandible (left) during biting. Practised users can dispense with this silicone register (for example, on the left) as the patient now already bites on the newly gained occlusion or bite elevation during biting. After fixating all occlusions, reocclusion is performed until the desired contact points have been attained.

DIRECT SYSTEM

OCCLUSIONVD

edelweiss
DENTISTRY